






One-Step Turn Lock



NEXXT MATRIXX® ALIF Integration



Multiple Insertion Options

PRODUCT FEATURES

- (A)** Integrated one-step turn lock feature to prevent Screw backout.
- (B)** Self-tapping, variable angle (0°-15° cephalad/caudal) Screws designed with tip-to-tail thread pattern for cancellous and cortical bone fixation.
- (C)** Anti-skid interface aids in Plate stabilization during placement and Screw fixation on vertebral body.
- (D)** Implant compatibility with **NEXXT MATRIX®** ALIF for single position procedural solution.
- (E)** Instrumentation designed to provide intra-operative flexibility with Plate and Construct Inserters.





Plate Insertion

Plate and Cage Construct Insertion


PLATE SPECS

Sizes




SCREW SPECS

Lengths (mm)




Diameters (mm)



CONSTRUCT SPECS


Clearance

Minimum of 2mm in between Ø5.5mm Screw



Angulation


Variable 0°-15°




INSTRUMENTATION

Self-Guided Solution


Intra-operative ease creating centered pilot holes



Locate



Dock



Plunge

Standard Instrument Set

Numerous handle options available upon request

Description
Hudson Axial Handle, Ratcheting
AO Axial Handle, Torque Limiting, 15inlb
ALIF Construct Inserter
ALIF Plate Inserter
Self-Guided Straight Awl
Self-Guided Angled Awl
Self-Guided Straight Drill
Straight Driver, T20 (2x)
Ball-Ended Lock Tool
Inserter Inner Shaft (2x)



Nexxt Spine, LLC
14425 Bergen Blvd, Suite B
Noblesville, IN 46060
(317)-436-7801
Info@NexxtSpine.com
NexxtSpine.com

For indications, contraindications, warnings, precautions, potential adverse effects and patient counselling information, see the package insert or contact your local representative; visit NexxtSpine.com for additional product information.

All rights reserved. All content herein is protected by copyright, trademarks and other intellectual property rights owned by Nexxt Spine, LLC and must not be redistributed, duplicated or disclosed, in whole or in part, without the expressed written consent of Nexxt Spine, LLC. This material is intended for health care professionals, the Nexxt Spine sales force and authorized representatives. Distribution to any other recipient is prohibited.