




InterMed Resources  
TENNESSEE

# Stone Collector

Improving stone fragment evacuation efficiency during RIRS.





# Stone Collector

1

## Ready-to-Use Configuration

The pre-installed tubing and integrated egress pressure-control knob allow for direct use.

2

## Multiple Capacity Options

Two collector volume configurations accommodate different clinical needs and procedural preferences.

3

## Bottle Handle for Secure Positioning

Allows the bottle to be held and positioned more easily during surgery.



4

## Reliable Specimen Transport


An additional sealed cap is provided to prevent leakage during specimen transport.



5


## Adjustable Negative Pressure Control

The negative aspiration pressure can be regulated using the pressure-control knob located on the egress tube.



### Separation and Filtration

Efficiently captures stone fragments  $\geq 1.2$  mm within the stone collector.



140 ml

180 ml

## Stone Collector Specification



Model	Volume	Item Code
180mL-L	180 ml	2521B-180L-PA001
140mL-S	140 ml	2521B-140S-PA002

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.