

When it comes to bipolar electrodes

# YOU HAVE A CHOICE.

High performance at a competitive price.

## **DORNIERDISC™** VAPORIZATION ELECTRODE

Dornier now offers a full line of validated bipolar electrodes compatible with market leading electrosurgical generators.\*

## RESECT & VAPORIZE

**DORNIER BIPOLAR ELECTRODES**  
Used in saline and designed for use in the prostate and bladder



\*Dornier bipolar saline electrodes are designed and validated to fit the Olympus bipolar WA22366A active and the WA22367A passive working elements. Use with Olympus UES 40 or ESG 400 generator systems.

<b>Product Name</b>	<b>Part Number</b>	<b>Product Description</b>
DornierDisc™	BIPDDSC	Dornier Plasma Disc Electrode
12 Degree Medium Cutting Loop	BIP12CLM	12 Degree Medium Cutting Loop Electrode
30 Degree Medium Cutting Loop	BIP30CLM	30 Degree Medium Cutting Loop Electrode
12 Degree Large Cutting Loop	BIP12CLL	12 Degree Large Cutting Loop Electrode
30 Degree Large Cutting Loop	BIP30CLL	30 Degree Large Cutting Loop Electrode
Needle Electrode	BIPNEEL	Needle Electrode
Medium Bladder Loop	BIPMBLA	Medium Bladder Loop Electrode

[www.dornier.com](http://www.dornier.com)

© 2019 DORNIER MEDTECH. ALL RIGHTS RESERVED. SUBJECT TO CHANGE WITHOUT NOTICE. DORNIERDISC™ IS A REGISTERED TRADEMARK OF DORNIER MEDTECH.

DMTA624-092019 REV A EN